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**Safety information for lithium batteries (CR2032/CR2430):**

1. Lithium batteries are lithium-manganese cells. If you are not intending to use the scales for a long period of time, you should remove the batteries from the device.
2. Never open the lithium batteries or throw them into fire or expose them to shocks, as harmful vapours could escape or batteries might explode.
3. Never touch leaking batteries with your bare hands.
4. If your eyes or hands are contaminated, rinse with a lot of water; in case of skin or eye irritation, consult a doctor.
5. Only dispose of cells that are completely discharged or packaged protected against polarity reversal in accordance with local disposal regulations.
6. Do not subject the cells to direct sunlight and extreme heat; risk of overheating!
7. Maintain a storage temperature of < 40 °C.


**Disposal of electrical and electronic devices**

The symbol on the product or its packaging indicates that this product must not be disposed of with regular household waste; instead, it must be disposed of at a collection point for the recycling of electrical and electronic equipment. For further details about recycling, contact your local council, municipal recycling companies or the shop where you bought the product.

Sincerely,

ADE (GmbH & Co.)

**Guarantee card – 3 years on the body analysis scales**

<ul style="list-style-type: none"> <li>• Sender</li> </ul> <hr/> <hr/> <hr/>	<ul style="list-style-type: none"> <li>• Reason for Claim</li> </ul> <hr/> <hr/> <hr/>
--	--

• Date of purchase

Please return this section and proof of purchase to your retailer.

**Model Jette / BA 1210  
Operating Instructions**

**Table of Contents**

1. How do the body analysis scales work? ..... Page 2
2. Conditions of Use ..... Page 2
3. First Use ..... Page 3
4. Weighing function (only weight measurement) ..... Page 4
5. Programming / Entering personal data ..... Page 4
6. Body analysis ..... Page 4
7. Error Messages ..... Page 6
8. Technical Data ..... Page 7
9. Warranty ..... Page 7

BODY FAT

BODY WATER

MUSCLE MASS

BONE MASS

ACTIVITY LEVEL

AUTO OFF

Thank you for purchasing these high-quality body analysis scales.

Please read these Operating Instructions carefully before using the product for the first time, and keep them in a safe place so that you can consult them again if you need further information.

## 1. How do the body analysis scales work?

With these new revolutionary scales you will receive information on your body composition regarding weight, percentage of body fat %, body water %, muscle mass %, and bone mass in kg. To start, you will need to enter your personal data one time in one of the memory locations: body height, age, gender, and your physical status (athletic/not athletic). Then weighing is as easy as on any normal scales. Then you simply press your individual memory button, step onto the scales, and in just a few seconds your current body status will be displayed. A small electrical current is used to measure the resistance. Water, which is mainly contained in muscles, is a good conductor, fat on the other hand blocks the current. Taking the measured resistance and the stored data into consideration, the body analysis scales can very precisely calculate your body composition.

## 2. Conditions of use

The measurements should take place without clothes, if possible, and always barefoot.

Since the body is subject to natural fluctuations (e.g. dehydration through sports or sauna, intake of food and fluids), constant conditions should be ensured when measuring, so as to receive comparable values.

Care should be taken that the measurements are always carried out at the same time of day. Directly after getting up in the morning, the water has not yet spread through the body. This means that the resistance is still relatively high, due to the bad conductivity in the body. With increasing movements throughout the day, the water spreads through the body, increasing conductivity, which causes the resistance, and thus the fat level too, to be shown as lower. To determine the actual changes of body composition, it is therefore important to measure always at same time. An easy way of doing this is by always measuring a quarter of an hour after getting up in the morning, either before or after using the toilet.

**Despite the minimum current, persons with cardiac pace makers or other implants, which emit frequencies, must not use the body analysis scales!** Other implants, such as a titanium plate in the knee, pose no danger. Metals, however, conduct electricity extremely well, which has the effect that body fat values are "beautified", meaning that the actual fat level is higher than shown by the scales. Tendencies, on the other hand, can be determined just as well.

## 8. Technical data

Load capacity x resolution:	Max. 180 kg x 100 g
Display resolution:	Percentage of body fat: 0.1 %
Display resolution:	Percentage of body water: 0.1 %
Display resolution:	Percentage of muscle mass: 0.1 %
Display resolution:	Bone mass 0.1 kg
Pers. memory locations:	8
Age:	10 – 100 years
Height:	100 – 220 cm
Dimensions:	318 mm x 294.5 mm x 27.5 mm
Batteries:	2 x CR 2032 (included)

Subject to technical changes as a result of product improvements.

## 9. Warranty

ADE agrees to rectify any material or manufacturing defects by offering free repair or replacement, for a period of 3 years from the purchase date. At the time of purchase, please get the retailer to fill in and stamp the guarantee card. In the case of guarantee claims, please return the scales to your retailer including guarantee card and stating the reasons for the claim.

**CE conformity. This device is interference-suppressed in accordance with the applicable EU Directive 2004/108/EC.**

Note: Extreme sources of electromagnetic radiation, such as those caused by the operation of radio equipment in the vicinity of the scales, can affect the value displayed. Once the interference has been eliminated, the product will once again function normally.

Note: the unit may need to be reset.

### Mandatory legal notice concerning battery disposal:

Do not dispose of batteries in domestic waste. As a consumer, you are legally required to return used batteries. You may return your used batteries to public collection points in your community or to any retailer that sells batteries of the same type.

### Note:

The following symbols are used on batteries:

Li = Battery contains lithium

Al = Battery contains alkaline

Mn = Battery contains manganese

CR (Li); AA (Al, Mn); AAA (Al, Mn)

Water makes up approx. 60 % of men's body mass. The value for women is about 55 % (due to the higher percentage of body fat). Of this total amount of body water, approx. two thirds are stored in the cells and is therefore called intracellular water. The remaining third is extracellular water.

Body fat, body water:

Women 					
Age	Body fat				Body water
	slim	normal	slightly overweight	overweight	normal
10-16	< 18 %	18-28 %	29-35 %	> 35 %	57-67 %
17-39	< 20 %	20-32 %	33-38 %	> 38 %	47-57 %
40-55	< 23 %	23-35 %	36-41 %	> 41 %	42-52 %
56-85	< 24 %	24-36 %	37-42 %	> 42 %	37-47 %
Men 					
Age	Body fat				Body water
	slim	normal	slightly overweight	overweight	normal
10-16	< 10 %	10-18 %	19-23 %	> 23 %	58-72 %
17-39	< 12 %	12-20 %	21-25 %	> 25 %	53-67 %
40-55	< 13 %	13-21 %	22-26 %	> 26 %	47-61 %
56-85	< 14 %	14-22 %	23-27 %	> 27 %	42-56 %

Bone mass:

	Women			Men		
Weight	< 45 kg	45-60 kg	> 60 kg	< 60 kg	60-75 kg	>75 kg
Bone	1.8 kg	2.2 kg	2.5 kg	2.5 kg	2.9 kg	3.2 kg

There are no general guide values for muscle mass.

## 7. Error messages

- "LO" = Battery low – Please replace the batteries with the correct battery type for your scales (2 x CR 2032).
- "----" = Overload – The scale's 180 kg maximum load capacity has been exceeded. Remove the object from the scales immediately to avoid damaging the scales' load cell unit.

In women, menstrual factors have an effect on the measurement results, since hormonal fluctuations influence the body's water balance. If the body contains all lot of water, the fat values are lower, if less water is contained, the values are a little higher.

When a woman is pregnant, there is no danger for the unborn child. But due to larger water deposits during pregnancy, the values are distorted making the measurement senseless.

ADE Germany does not assume liability for damages or losses, caused by the body analysis scales, or for claims of third parties. This product is solely intended for private use. It is not designed for professional use in hospitals or medical facilities.

## 3. First use

**Safety guidelines:** This device is not to be used by persons (including children) with limited physical, sensory or mental abilities or who are lacking experience and understanding, unless they are supervised by persons, who are responsible for their safety or have been instructed in the use of the device. It should be ensured that children do not play with the device.

Not suitable for persons with electronic implants (cardiac pace makers, etc.).

Place the scales on a flat and solid surface (tiles, parquet, etc.). Incorrect measurements make occur on carpets. Stand in the centre of the scale surface to prevent tipping.

Attention! There is a danger of slipping if the scales' surface is wet or if feet are wet.

Ensure that no objects fall onto the scales; otherwise, the glass might break.

Use only a damp cloth to clean the scales. Never use solvents or abrasive cleaning agents. Do not immerse in water.

**Preparation:** Open the battery compartment cover on the underside of the scales and remove the battery safety strip. Please check for correct polarity of the batteries (+/-).

On the underside of the scales, a small push button "UNIT" is provided, which will allow country specific settings (kg = kilogram / lb = pound). Please ensure that the scales are configured to show the unit of weight that you normally use.

#### 4. Weighing function (only weight measuring)

Simply step onto the scales. Stand still, do not hold on to anything. Your weight will be displayed as soon as the scales have stabilised.

#### 5. Programming / Entering personal data

Enter and store your personal data (female or male, athletic or non-athletic, body height, age) in one of the memory locations P1 – P8 (for up to 8 individuals) before using the scales for the first time. The input buttons are located below the display. To enter your personal data, touch the buttons: SET (memory button), ◀ (minus button), ▶ (plus button), and complete the individual steps (quick number cycle when touching and holding the button for more than 1 second):

1. Touch the SET button. The memory location selection display (e.g. P1) will flash. Touch the ◀ or ▶ button to select the desired memory location (P1 to P8), and save by touching the SET button.
2. After you have selected the memory location (e.g. memory location P3), the gender icon flashes. Press the ◀ or ▶ button to select the correct gender (♂ male, ♀ female, male athletic ♂ or ♀ or female athletic ♀), then save by pressing the SET button.
3. Once you have selected your gender and athletic status, the body height display will begin to flash (default setting 160 cm). Touch the ◀ or ▶ button to enter your body height, and save by pressing the SET button.
4. After you have entered your body height (e.g. 166 cm), the age display will start flashing (default setting 30). Touch the ◀ or ▶ button to enter your age, and save by pressing the SET button.

#### 6. Body analysis

An analysis is only possible with bare feet. The personal data must be entered beforehand (see page 4).

1. For your first analysis, select your personal memory location via the SET button. Carefully place your bare feet on the stainless steel strips on the scales as soon as 0.0 is displayed, and stand still.

2. For any subsequent measurement, simply step onto the scales, carefully placing your bare feet on the stainless steel strips of the scales. The scales will activate automatically.
3. Your weight will be displayed first, and the scales will then search for your personal memory location. Should two memory locations with very similar weight values exist, then both memory places (e.g. P1 and P2) will be displayed simultaneously. Select the memory location on the left by pressing the ◀ button, or the one on the right by pressing the ▶ button. The signal 'o' will now run twice from left to right across the display, while the body analysis is calculated. Keep standing still. Once complete, the result values will be shown on the display as follows:

Constant Your weight in kg

1. Your body fat in % (FAT)\*\*\*
2. Your body water in % (TBW)
3. Your muscle mass in % (MUS)
4. Your bone mass in kg (BONE)

\*\*\* An analysis will be shown at the very bottom of the display, while your body fat percentage is displayed:

	= Body fat % too low
	= Body fat % in healthy range
	= Increased body fat %
	= Body fat % too high

The measured values will be displayed twice more in succession, then the scales will deactivate automatically.

The percentage of body fat shows the percentage of fat in relation to the total weight of a body. There is no general recommendation for this value, since this depends significantly on gender (due to the different physique) and age. Body fat is vital for the daily body functions. It protects organs, cushions joints, regulates body temperature, stores vitamins and acts as an energy storage for the body. With the exception of high-performance athletes, a body fat percentage for women between 10 % and 13 % and for men between 2 % and 5 % is regarded as vital.

The percentage of body fat is not a clear indicator for health. A body weight and a percentage of body fat which are too high are related to many diseases of civilisation such as diabetes, cardiovascular diseases, etc. A lack of exercise and the wrong diet are the most common causes for these diseases. There is a clear causal connection between being overweight and a lack of exercise.

A significant part of the human body consists of water. This body water has a different distribution throughout the body. Lean muscle tissue contains about 75 % water, blood approx. 83 % water, body fat approx. 25 % water and the bones have a water percentage of about 22 %.

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