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Statutory notice of liability for disposal of batteries

Batteries are not domestic refuse. As a consumer, you are legally liable to return used batteries. You can return your used batteries to the public collection points in your area or anywhere where batteries of the respective type are sold.

Note:

You will find these markings on the batteries:

- Li = The battery contains lithium
- Al = The battery contains alkali
- Mn = The battery contains manganese


CR (Li); AA (Al, Mn); AAA (Al, Mn)

Disposing of used electrical and electronic equipment

The symbol on the product or its packaging indicates that this product must not be treated as normal household refuse; rather, you must dispose of it at a collection point for recycling electrical and electronic equipment. You can get further information from your local council, the municipal recycling authorities or the shop where you bought the product.

Yours truly
ADE GmbH & Co.

Guarantee card – 3 years on the body analysing scale

<ul style="list-style-type: none"> • Sender <hr/> <hr/> <hr/>	<ul style="list-style-type: none"> • Reason for complaint <hr/> <hr/> <hr/> <hr/>
--	--

• Date of purchase

Please return this guarantee card together with the receipt to your retailer.

**Model Tanja / BA 719
Operating Instructions**

Table of Contents

- 1. How does the body analyser work?..... page 2
- 2. Conditions to use..... page 2
- 3. Commissioning..... page 3
- 4. Weighing function (weight measurement only)..... page 4
- 5. Programming / enter of personal data..... page 4
- 6. Body Analysis..... page 5
- 7. Error Messages..... page 7
- 8. Technical Data..... page 7
- 9. Guarantee..... page 7

BODY FAT

BODY WATER

MUSCLE MASS

STEP ON
FUNCTION

AUTO OFF

Thank you for purchasing this high-quality product made by ADE. You have acquired a set of scales that meet the highest technical standards and are of proven quality.

Please read these operating instructions carefully before using the product for the first time and keep them in a safe place where you can find them whenever you need the information contained in the instructions.

1. How does the body analyser work?

Through this new revolutionary scale, you will be able to get information about your body composition with reference to your weight, body fat %, body water % and muscle mass %. First, you have to enter once your personal data like body height, age, gender under your selected personal memory location. Afterwards a measurement is as simple as on a normal bathroom scale. You have only to push your personal memory location, step on the scale and within a few seconds, your measured values will be indicated. The analysis based on the measurement of the body's electrical resistance. A low, safe electrical signal passes freely through fluids contained in lean tissue, such as muscle and blood, but meets resistance passing through fat tissue. The body analyser accurately measures the resistance and with reference to your personal data use it to calculate elements of body composition.

2. How to use

If possible, the measurement should be done undressed but always barefoot.

By nature, a body is subject of fluctuations (for example loss of water after taking exercise or after a sauna visit, after a meal, etc.). To enable comparable readings it is important to take the measurements as possible always under the same conditions.

In order to ensure that the results of analysis are accurate and consistent, please keep the measurement conditions constant always at the same daytime. In the morning, directly after you've woken up, your body will be dehydrated and also the majority of fluid will be stored in the central trunk area. This means that the resistance is relatively high because of the poor conductivity in the body. By intensify moving during the day progresses the fluid becomes more evenly distributed through the limbs and increasing the accuracy of the readings. To get a realistic result it is important to take your readings at the same time and under the same conditions every day. The easiest way to stick the rule, measure always in the morning, 15 minutes after you've got up, either before or after you've gone to the bathroom.

In spite of a very low electrical signal, this body analyzer is not suitable for persons with heart pacemaker or other electronic medical implants! There is no risk for persons with none electronic implant for instance a titanium plate in the knee. However, metal conducts the electric current extremely well and therefore the scale implicates the body fat value "beautified" – means the actual body fat value is higher as indicated. However, tendencies may be determined as well.

7. Error Messages

"FFFF"	=	Overload display – The 150 kg maximum capacity of the scales has been exceeded. Immediately remove the load from the scales; otherwise, the load cell could be damaged.
"Lo"	=	Indication that the battery is flat – Please replace the batteries with four 1.5 V type AA.

8. Technical Data

Capacity x graduation:	max. 150 kg x 100 g
Display resolution:	body fat: 0.1%
Display resolution:	body water: 0.1%
Display resolution:	muscle mass: 0.1%
Personal memory loc.:	12 positions
Age:	9 – 99 years
Body height:	90 – 250 cm
Dimensions:	330 mm x 330 mm x 36 mm
Batteries:	4 x 1.5V AA (included)

Subject to technical changes due to further-development of the product.

5. Guarantee

ADE guarantees to rectify any defects due to material or manufacturing defects by means of repair or replacement free of charge for a period of three years. At the time of purchase, please get the retailer to fill in and stamp the guarantee card. In the case of guarantee claims, please return the scales to your retailer including the guarantee card and stating the reasons for the claim.

CE conformity

This device is interference-suppressed in accordance with the applicable EU directive 2004/108/EC

Note: Extreme electromagnetic influences, e.g. when operating radio communication equipment is being operated next to the scales, can affect the displayed value. When the interference is over, you can use the product for its intended purpose again; it may be necessary to switch it on again.

The body fat content shows the amount in percentage of stored fat from the total weight of a body. There is no generally accepted set value as this varies significantly according to gender (due to the different build) and age. Body fat is vital for basic body functions. It's protecting organs, cushioning joints, regulating body temperature, storing vitamins and serves the body as energy storage.

As essential for women a minimum body fat proportion between 10% and 13% and for men between 10% and 13% is requested - with the exception of competitive athletes.

Body fat is not a clear indicator for health. Too high body weight and a too high body fat proportion are linked to the incidence of many civilization diseases such as diabetes, cardiovascular diseases, etc. Lack of movement and improper nutrition are often the prime course for these diseases. There is a obvious direct relationship between overweight and lack of movement.

A significant part of the human body consists of body water. This body water is differently distributed throughout the whole body. Lean muscle tissue contains approx. 75% body water, blood approx. 83% body water, body fat approx. 25% body water and bones have a water content of about 22%.

The average total body water percentage range for healthy men is between 50 – 65%. The equivalent value for women is approximately 45 – 60% due to the higher body fat content. About two third of the total body water is stored in the cells and is therefore called intra cellular. The other third of the total body water is extra cellular body water.

Women 					
AGE	Body Fat				Body Water
	slim	normal	slightly overweight	overweight	normal
10-16	<18%	18-28%	29-35%	>35%	57-67%
17-39	<20%	20-32%	33-38%	>38%	47-57%
40-55	<23%	23-35%	36-41%	>41%	42-52%
56-85	<24%	24-36%	37-42%	>42%	37-47%

Men 					
AGE	Body Fat				Body Water
	slim	normal	slightly overweight	overweight	normal
10-16	<10%	10-18%	19-23%	>23%	58-72%
17-39	<12%	12-20%	21-25%	>25%	53-67%
40-55	<13%	13-21%	22-26%	>26%	47-61%
56-85	<14%	14-22%	23-27%	>27%	42-56%

There are no validated standards for muscle mass.

Due to menstrual cycle and hormonal changes women experience greater changes in their hydration level. Retaining fluid may cause the weight to fluctuate day-to-day. This will also affect the body fat percentage readings. If the hydration level is high, the body fat percentage is slightly less – a less hydration level means the body fat percentage is slightly higher.

There is no healthy risk for pregnancy. The body analyzer is absolutely safe for the mother and her unborn baby. However, due to the changes in water content in the mother's body, the readings should not be interpreted as completely accurate.

ADE Germany takes no responsibility for any damages or losses caused by the use of this body analyzer scale or for Third Party claims. This product is exclusively intended for private use. This product is not designed for commercial or professional use in hospitals or other medical facilities.

3. Commissioning

Safety guidelines:

Not suitable for persons with electronic implants (heart pacemaker, etc.).

Place the scale on a plane and firm surface (tiles, hardwood flooring, etc.) to avoid a wrong measurement on carpeted floor.

Always stay in the centre of the scale to prevent the scales tipping to the side.

Attention - Hazard of slips! Do not stand on the scale with wet feet or on a slippery surface of the weighing platform.

Ensure that no objects fall onto the scale; the glass surface is breakable.


Clean the scale only by using a slightly damp cloth. Never use solvents or abrasive cleaners. Do not submerge the scale under water.

Preparation:

Open the lid of the battery compartment on the bottom of the scales and insert the four supplied 1.5 V AA batteries. Please observe the correct polarity of the batteries (+/-).

4. Weighing function (weight measurement only)


1. Switch on the scale by using the step-on-function. Just step on the scale. Remain stationary on the weighing platform and do not hold you tight. After a few seconds, your weight will be indicated.
2. After weighing, step off the scale. The scale will switch off automatically after a few seconds.


5. Programming / Enter of personal data

Before first time use your personal data (gender: female or male, body height, age) must be programmed and saved under selected locations P-01 up to P-12 (for up to 12 persons). The enter keys are located at the head front of the scale. To enter your personal data please utilise the buttons: ◀ (minus button), ▶ (plus button), SET (save button) and follow the particular points:


1. Place the scale on a table. Switch on the scale by pressing the SET-button, the display indicates a flashing "P-01". To select your personal location (P-01 – P-12), press the ◀ or ▶ button. After the selection of your personal location, press the SET-button to save the reading.


2. After you selected your personal location, the symbol for gender starts to flash. To select the requested gender press the ◀ or ▶ button. After the selection of your gender press the SET-button to save the reading.


3. After you selected your gender, the symbol for body height starts to flash. To enter your heights press the ◀ or ▶ button. After the selection for your height press the SET-button to save the reading.

4. After you selected your body height, the symbol for age starts to flash. To enter your age press the ◀ or ▶ button. After the selection of your age press the SET-button to save the reading.

5. After input completion, the display indicates again your entered location (the entered personal location can be changed by using the ◀ or ▶ button) and your entered personal data. A few seconds later, the scale automatically powers off. To enter personal data of further persons please repeat the steps listed below.


6. Body Analysis

An analysis is only possible in bare feet. Your personal data must be entered first (please see page 4 - 5).

1. Press the SET-button and select your personal location by using the ◀ or ▶ button.
2. At first, your personal data will be displayed (age, body height, and gender). Shortly after that, „0.0“ will be indicated. As soon as the display shows „0.0“ the scale is prepared for the measurement.
3. Carefully step on the scale in bare feet and remain stationary while the measurement is now taking place.
4. After completion of the measurement will appear in a sequence:
 1. your body fat, your actual weight, your body height, and your age
 2. your body water, your ideal weight, and your actual calories
 3. your muscle mass, your bone mass, and your daily-required calories

EXAMPLE:

5. After weighing, step off the scale. The scale switches off automatically after a few seconds.

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