**Hinweis:**

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**Note:**

You will find these markings on the batteries:

Li = The battery contains lithium

Al = The battery contains alkali

Mn = The battery contains manganese

CR (Li); AA (Al, Mn); AAA (Al, Mn)


**Disposal of old electrical and electronic equipment**

The symbol on the product or its packaging indicates that this product must not be treated as normal household refuse; rather, you must dispose of it at a collection point for recycling electrical and electronic equipment. You can get further information from your local council, the municipal recycling authorities or the shop where you bought the product.

Yours truly

ADE GmbH & Co.

**Guarantee card – 3 years guarantee for the body analysing scale**

• Sender

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

• Reason for complaint

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

• Date of purchase

Please return this guarantee card together with the receipt to your retailer.

BODY FAT

BODY WATER

MUSCLE MASS

BONE MASS

CALORIE REQUIREMENT

**Model Marike / BA 1100  
 Instruction Manual**

**Directory**

1. How does the body analyser work?.....Page 2  
 2. Conditions to use.....Page 2  
 3. Commissioning.....Page 3  
 4. Programming / enter of personal data.....Page 4  
 5. Body Analysis.....Page 4  
 6. Automatic Person Recognition.....Page 5  
 7. Error Messages.....Page 7  
 8. Technical Data.....Page 7  
 9. Warranty.....Page 7

Thank you for purchasing this high-quality product made by ADE. You have acquired a body analyser scale that meet the highest technical standards and are of proven quality. Please read these operating instructions carefully before using the product for the first time and keep them in a safe place where you can find them whenever you need the information contained in the instructions.

## 1. How does the body analyser work?

Through this new revolutionary scale, you will be able to get information about your body composition with reference to your weight, body fat %, body water %, muscle mass %, bone mass kg, daily calorie requirement (basal metabolic rate) in kcal. First of all, you have to enter once your personal data like gender (male, female, or male athlete, female athlete), age and body height under your selected personal memory location. Afterwards a measurement is as simple as on a normal bathroom scale. Afterwards a measurement is as simple as on a normal bathroom scale. You have only to select your personal memory location, step on the scale and within a few seconds, your measured values will be indicated. The analysis based on the measurement of the body's electrical resistance. A low, safe electrical signal passes freely through fluids contained in lean tissue, such as muscle and blood, but meets resistance passing through fat tissue. The body analyser accurately measures the resistance and with reference to your personal data use it to calculate elements of body composition.

## 2. How to use

If possible, the measurement should be done undressed but always barefoot.

By nature, a body is subject of fluctuations (for example, loss of water after taking exercise or after a sauna visit, after a meal, etc.). To enable comparable readings it is important to take the measurements as possible always under the same conditions.

In order to ensure that the results of analysis are accurate and consistent, please keep the measurement conditions constant always at the same daytime. In the morning, directly after you've woken up, your body will be dehydrated and also the majority of fluid will be stored in the central trunk area. This means that the resistance is relatively high because of the poor conductivity in the body. By intensify moving during the day progresses the fluid becomes more evenly distributed through the limbs and increasing the accuracy of the readings. To get a realistic result it is important to take your readings at the same time and under the same conditions every day. The easiest way to stick the rule, measure always in the morning, 15 minutes after you've got up, either before or after you've gone to the bathroom.

**In spite of a very low electrical signal, this body analyser is not suitable for persons with heart pacemaker or other electronic medical implants!** There is no risk for persons with none electronic implant for instance a titanium plate in the knee. However, metal conducts the electric current extremely well and therefore the scale implicates the body fat value "beautified" – means the actual body fat value is higher as indicated. However, tendencies may be determined as well.

## 7. Error Messages

„Err“ = Overload display - The 180 kg maximum capacity of the scale has been exceeded. Please immediately remove the critical load from the scale; otherwise the load cell could be damaged.

„Lo“ = Indication the batteries are flat - Please replace the batteries with the intended type (4 x 1.5V AA8).

## 8. Technical Data

Capacity x graduation:	max. 180 kg x 100 g
Display resolution:	body fat: 0.1%
Display resolution:	body water: 0.1%
Display resolution:	muscle mass: 0.1%
Display resolution:	bone mass: 100 g
Personal memory loc.:	10 positions
Age:	8 – 100 years
Height:	80 – 240 cm
Dimensions:	approx. 361 mm x 241 mm x 23 mm
Batteries:	four x 1.5V AAA (included)

Subject to technical changes due to further-development of the product.

## 9. Guarantee

ADE guarantees to rectify any defects due to material or manufacturing defects by means of repair or replacement free of charge for a period of three years. In the case of guarantee claims, please return the scale together with the guarantee card to your retailer and stating the reasons for the claim.

**CE conformity**  
**This device is interference-suppressed in accordance with the applicable EU directive 2004/108/EC**

Note: Extreme electromagnetic influences, e.g. when operating radio communication equipment is being operated next to the scales, can affect the displayed value. When the interference is over, you can use the product for its intended purpose again; it may be necessary to switch it on again.

### Statutory notice of liability for disposal of batteries:

Batteries are not domestic refuse. As a consumer, you are legally liable to return used batteries. You can return your used batteries to the public collection points in your area or anywhere where batteries of the respective type are sold.

Body fat is not a clear indicator for health. Too high body weight and a too high body fat proportion are linked to the incidence of many civilization diseases such as diabetes, cardiovascular diseases, etc. Lack of movement and improper nutrition are often the prime course for these diseases. There is an obvious direct relationship between overweight and lack of movement.

A significant part of the human body consists of body water. This body water is differently distributed throughout the whole body. Lean muscle tissue contains approx. 75% body water, blood approx. 83% body water, body fat approx. 25% body water and bones have a water content of about 22%.

The average total body water percentage range for healthy men is between 50 – 65%. The equivalent value for women is approximately 45 – 60% due to the higher body fat content. About two third of the total body water is stored in the cells and is therefore called intra cellular. The other third of the total body water is extra cellular body water.

For men is the normal range of body mass index between 20 and 25, for women 19 to 24. When people in the second half of life a body mass index to 26 is considered as still normal.

Women 					
AGE	Body Fat				Body Water
	slim	normal	slightly overweight	overweight	normal
10-16	<18%	18-28%	29-35%	>35%	57-67%
17-39	<20%	20-32%	33-38%	>38%	47-57%
40-55	<23%	23-35%	36-41%	>41%	42-52%
56-85	<24%	24-36%	37-42%	>42%	37-47%
Men 					
AGE	Body Fat				Body Water
	slim	normal	slightly overweight	overweight	normal
10-16	<10%	10-18%	19-23%	>23%	58-72%
17-39	<12%	12-20%	21-25%	>25%	53-67%
40-55	<13%	13-21%	22-26%	>26%	47-61%
56-85	<14%	14-22%	23-27%	>27%	42-56%

There are no validated standards for muscle mass.

Due to menstrual cycle and hormonal changes women experience greater changes in their hydration level. Retaining fluid may cause the weight to fluctuate day-to-day. This will also affect the body fat percentage readings. If the hydration level is high, the body fat percentage is slightly less – a less hydration level means the body fat percentage is slightly higher. There is no healthy risk for pregnancy. The body analyzer is absolute safe for the mother and her unborn baby. However, due to the changes in water content in the mother's body, the readings should not be interpreted as completely accurate.

ADE Germany takes no responsibility for any damages or losses caused by the use of this body analyzer scale or for Third Party claims. This product is exclusive intended for private use. This product is not designed for commercial or professional use in hospitals or other medical facilities.

### 3. Commissioning

**Safety guidelines:** The appliance is not intended for use by persons (including children) with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, unless a person responsible for their safety has given them supervision or instruction concerning use of the appliance. Children should be supervised to ensure that they do not play with the appliance.

Not suitable for persons with electronic implants (heart pacemaker, etc.).

Place the scale on a plane and firm surface (tiles, hardwood flooring, etc.) to avoid a wrong measurement on carpeted floor

Attention - Hazard of slips! Do not stand on the scale with wet feet or on a slippery surface of the weighing platform.

Please care that no heavy object may fall on the scale. The glass surface is breakable.

Clean the scale only by using a slightly damp cloth. Do not use solvents or abrasive cleaning agents. Do not submerge the scale under water.

#### Preparation:

Open the lid of the battery compartment on the bottom of the scales and remove the battery safety strip. Please observe the correct polarity of the batteries (+/-).

To change the country specific setting (kg = kilograms / lb = pound or st: lb = stone: pound) switch-on the scale by touching the SET-button and select the requested country specific unit of weight by touching the ◀-button. Please pay attention to use the scale under consideration of your country specific unit of weight.

#### 4. Programming / Enter of personal data

Before first time use, your personal data like gender (female or male), athlete, age, and body height must be programmed and saved under selected locations P1 – P9 + P0 guest (up to 9 persons + 1 guest). The scale is equipped with highly sensitive touch buttons. The touch buttons are located below the display. To enter your personal data please utilise the buttons: SET (on and save button), ◀ (minus button), ▶ (plus button) and follow the particular points:

1. Lightly touch the SET-button two times. The display indicates flashing “P01. To select your personal location (P0 – P9) touch the ◀- or ▶-button. After the selection of your personal location, touch the SET-button to save the reading.
2. After you have selected your personal location (e.g. P3), the symbol for gender starts to flash. To select the requested gender (♂ male, ♀ female, or ♂ female athlete, ♀ male athlete) touch the ◀- or ▶-button. After the selection of your gender, touch the SET-button to save the reading.
3. After you have selected your gender (e.g. ♀ male), the symbol for age starts to flash (default age 30). To enter your age, touch the ◀- or ▶-button. After the selection of your age, touch the SET-button to save the reading.
4. After you have selected your age, the symbol for height starts to flash (default 170 cm). To enter your height, touch the ◀- or ▶-button. After the selection of your height, touch the SET-button to save the reading.

After input completion, the display indicates your selected personal data and “0.0 kg”. The scale is ready to use and you may now start your body analysis. Without further activity the scale automatically power-off. Your personal data are stored.

If you push the ◀- or ▶-button within 2 seconds after the first complete input, you can directly enter the personal data of a further person. Alternatively, you repeat the complete steps listed below for further entries of personal data.

#### 5. Body Analysis

An analysis is only possible in bare feet. Your personal data must be entered first (please see page 4).

1. Switch-on the scale for an initial analysis by touching the SET-button and select your personal location by using the ◀- or ▶-button. The display indicates your personal data and “0.0 kg”.

2. Now, carefully step in bare feet on the right and left stainless steel electrodes (measuring units) of the scale and remain stationary.

3. At first, your weight will be indicated twice on the display. Remain stationary on the scale. Shortly after, while the analysis is taking place, the signal “o” starts to move on the display. After completion of the measurement, your measured values will be shown two times successively in each of 2 ads for 7 seconds as follows:

1 <sup>st</sup> . ad	your weight in kg your body fat in % (Fat) your body water in % (Water)	
2 <sup>nd</sup> . ad	your daily calorie requirement in kcal (BMR) your bone mass in kg (Bone) your muscle mass in % (Muscle)	

Afterwards the scale switches-off automatically.

#### 6. Automatic Person Recognition

After the first measurement the scale have a person recognition if the balance has already saved once your values and your weight didn't changed more than ± 2 kg since the last weighing.

1. Simply, carefully step in bare feet on the right and left stainless steel electrodes (measuring units) of the scale and remain stationary. The scale will automatically switch-on and only your weight will be displayed
2. Shortly after, while the analysis is taking place, the signal “o” starts to move on the display. Remain stationary on the scale. After completion of the measurement, your measured values will be shown two times successively in each of 2 ads for 7 seconds as described above.
3. Should be saved multiple users with similar values, the scale displays the two most similar memory locations for selection.
4. If the scale cannot allocate the measured values to any location, these are stored on the memory P0.

The body fat content is shown as the amount of stored fat as a percentage of total body weight. There is no generally accepted set value as this varies significantly according to gender (due to the different build) and age. Body fat is vital for basic body functions. It's protecting organs, cushioning joints, regulating body temperature, storing vitamins and serves the body as energy storage. As essential for women a minimum body fat proportion between 10% and 13% and for men between 10% and 13% is requested - with the exception of competitive athletes.

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